


For further assistance on how to
use your Tenovi device,
please contact a Tenovi representative
at:

Phone: (714) 418-5658

**Email:
support@tenovi.com**



- 1** Plug your Gateway into a power outlet with the provided power cord and wait until the LED on the Gateway turns red.



- 2** Take a Daily Tenovi Device Measurement

- Please read the instruction manual of each device prior to taking a measurement.
- Take a measurement using the Tenovi Device(s) included with the Gateway.
- The Gateway LED will turn yellow while a measurement is being taken and will flash green when the measurement has been successfully transmitted.
- If all measurements prescribed by your doctor have been taken, the LED on the Gateway will change to Green for the remainder of the day.

3 If the LED ring turns yellow, this indicates that the measurement has not been properly transmitted. If it remains yellow for more than 15 minutes, try the gateway to another location with a better signal. If the issue persists, please call our customer service team.



If the Gateway flashes blue for more than 30 seconds, wait 5 seconds and plug the power adapter back into the Gateway. If the issue persists, please get in touch with customer service.



5 If the LED ring is Red, you must take a measurement for the day. Unless otherwise recommended by your physician, you should take one measurement per day.



Measurement Needed

4 If the LED ring is Green, you have properly taken a measurement for the day, and that measurement has been successfully transmitted.



6 If the LED ring flashes Red after you take a reading, you have triggered a physician-assigned alert. If this happens, take one more reading. If the LED ring flashes RED again, we recommend you write down the measurement value, and call your doctor, and inform them that you have triggered an alert with the value recorded.



Flashing Red